

UROFAST URINE REAGENT STRIPS



ONE-PARAMETER

Blood pH
Glucose Ketone
Protein Ascorbic Acid

TWO-PARAMETERS

Glucose/Protein Nitrite/Leukocyte
Glucose/Ketone Microalbumin/Creatinine
pH/Specific Gravity Glucose/Protein

THREE-PARAMETERS

Glucose/Protein/pH

FOUR-PARAMETERS

Glucose/Protein/Ph/Blood
Glucose/Protein/pH/Specific Gravity
Leukocyte/Nitrite/Protein/Glucose

FIVE-PARAMETERS

Blood/Ketone/Glucose/Protein/pH
Leukocyte/Nitrite/Blood/Protein/Glucose

SIX-PARAMETERS

Protein/pH/Blood/Ketone/Bilirubin/Glucose

SEVEN-PARAMETERS

Urobilinogen/Protein/pH/Blood/Ketone/
Bilirubin/Glucose

EIGHT-PARAMETERS

Glucose/Protein/pH/Ketone/Urobilinogen/
Bilirubin/Blood/Nitrite

NINE-PARAMETERS

Glucose/Protein/pH/Ketone/Urobilinogen
Bilirubin/Blood/Nitrite/Specific Gravity

TEN PARAMETERS

Glucose/Protein/pH/Bilirubin/Blood/Ketone Urobilinogen/Ketone/Nitrite/
Specific Gravity/Leukocyte

ELEVEN PARAMETERS

Glucose/Protein/pH/Bilirubin/Blood/Ketone Urobilinogen/Ketone/Nitrite/
Specific Gravity/Leukocyte/Ascorbic Acid



UROFAST-101 (ideal for physician's office or small lab) is a semi-automated and CLIA-waived 60 tests per hour urine strip reader with 1000 test memory, space-saving design has patient ID, and LIS capability, flip screen, and 110/220 volt design.

UROFAST-720 (ideal for small lab) is a semi-automated and CLIA-waived 500 tests per hour urine strip reader with 2000 test memory, space-saving design has patient ID, and LIS capability, flip screen, and 110/220 volt compatibility